

Neuropeptide Y Antibody FITC Conjugated

Catalog No: #C00053F

Package Size: #C00053F 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuropeptide Y Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 36-48 97 derived from human NPY
Conjugates	FITC
Target Name	Neuropeptide Y
Other Names	PYY4; Pro-neuropeptide Y; NPY
Accession No.	Swiss-Prot#P01303NCBI Gene ID4852
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.