

Caspase-6 (NT) Antibody Biotin Conjugated

Catalog No: #C00063B



Package Size: #C00063B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-6 (NT) Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 24-56 293 derived from human Caspase-6 subunit p18
Conjugates	Biotin
Target Name	Caspase-6 NT
Other Names	MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2; CASP6
Accession No.	Swiss-Prot#P55212NCBI Gene ID839
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves poly(ADP-ribose) polymerase in vitro, as well as lamins. Overexpression promotes programmed cell death.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.