

# Caspase-6 (NT) Antibody HRP Conjugated

Catalog No: #C00063H

Package Size: #C00063H 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Caspase-6 (NT) Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 24-56 293 derived from human Caspase-6 subunit p18
Conjugates	HRP
Target Name	Caspase-6 NT
Other Names	MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2; CASP6
Accession No.	Swiss-Prot#P55212NCBI Gene ID839
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

## Background

Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves poly(ADP-ribose) polymerase in vitro, as well as lamins. Overexpression promotes programmed cell death.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.