

## CCBR1 Antibody HRP Conjugated

Catalog No: #C06665H



Package Size: #C06665H 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CCBR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 211-234 2501 derived from human xCT CCBR1
Conjugates	HRP
Target Name	CCBR1
Other Names	xCT; CCBR1; Cystine glutamate transporter; Amino acid transport system xc-; Calcium channel blocker resistance protein CCBR1; Solute carrier family 7 member 11; SLC7A11
Accession No.	Swiss-Prot#Q9UPY5NCBI Gene ID23657
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

## Background

Sodium-independent, high-affinity exchange of anionic amino acids with high specificity for anionic form of cystine and glutamate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.