

CAPD3 Antibody HRP Conjugated

Catalog No: #C06981H

Package Size: #C06981H 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CAPD3 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 670-715 1498 derived from human CAPD3
Conjugates	HRP
Target Name	CAPD3
Other Names	CAPD3; hcp-6; CAP-D3; hHCP-6; hCAP-D3; Condensin-2 complex subunit D3; Non-SMC condensin II complex subunit D3; NCAPD3; KIAA0056
Accession No.	Swiss-Prot#P42695NCBI Gene ID23310
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Regulatory subunit of the condensin-2 complex, a complex which establishes mitotic chromosome architecture and is involved in physical rigidity of the chromatid axis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.