

GPCPD1 Conjugated Antibody

Catalog No: #C31870

Package Size: #C31870-Conjugated 50ul #C31870-AF680 100ul #C31870-Biotin 100ul #C31870-AF555 100ul #C31870-AF488 100ul
#C31870-AF594 100ul #C31870-AF647 100ul #C31870-AF750 100ul #C31870-AF350 100ul #C31870-AF405 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPCPD1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Fusion protein of human GPCPD1
Conjugates	AF680, AF594, AF647, AF750, AF350, AF405
Target Name	GPCPD1
Other Names	EDI3; GDE5; GDPD6; PREI4
Accession No.	Swiss-Prot#: Q9NPB8NCBI Protein#: BC027588
Calculated MW	76 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 500-2000

IF: 1:50-1:200

Background

May be involved in the negative regulation of skeletal muscle differentiation, independently of its glycerophosphocholine phosphodiesterase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.