

β-Actin Rabbit Conjugated Antibody

Catalog No: #C52901



Package Size: #C52901-AF350 100ul #C52901-AF405 100ul #C52901-AF488 100ul #C52901-AF555 100ul #C52901-AF594 100ul #C52901-AF647 100ul #C52901-AF680 100ul #C52901-AF750 100ul #C52901-Biotin 100ul #C52901-Conjugated 50ul

Orders: order@signalwayantibody.com  
Support: tech@signalwayantibody.com

Description

Product Name	β-Actin Rabbit Conjugated Antibody
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	beta Actin Antibody detects endogenous levels of total beta-Actin protein.
Immunogen Description	A synthesized peptide derived from human beta Actin
Conjugates	AFθ• §ε η´
Other Names	ACTIN; ACTB; BETA-ACTIN; beta actin; actin beta;
Calculated MW	42kDa
Formulation	Liquid
Storage	Store at 4Λ C in dark for 6 months

Application Details

Suggested Dilution:
AF350 conjugated: most applications: 1: 50 - 1: 250
AF405 conjugated: most applications: 1: 50 - 1: 250
AF488 conjugated: most applications: 1: 50 - 1: 250
AF555 conjugated: most applications: 1: 50 - 1: 250
AF594 conjugated: most applications: 1: 50 - 1: 250
AF647 conjugated: most applications: 1: 50 - 1: 250
AF680 conjugated: most applications: 1: 50 - 1: 250
AF750 conjugated: most applications: 1: 50 - 1: 250
Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Published Papers

Wang Hechen, Shen Xudan, Liu Jiatong, Zhu Xinlan, Zeng Su, Cai Sheng el at., Enhancing HepG2 cell apoptosis with a combined nanoparticle delivery of miR-128β 3p agomir and Oroxin B: A novel drug delivery approach based on PI3KB β ¶, Asian journal of pharmaceutical sciences, (2025)

PMID:40487120

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.