

VEGFR2 Polyclonal Antibody

Catalog No: #C90236



Package Size: #C90236-1 50ul #C90236-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VEGFR2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;FC
Species Reactivity	Human;Mouse;Rat
Conjugates	Unconjugated
Calculated MW	151kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P□1:100-500

Flow-Cyt□1&mu;g/Test

Published Papers

Esther J. J. Habets;Eef J. Hendriks;Martin J. B. Taphoorn;Linda Douw;Aeilko H. Zwinderman;W. Peter Vandertop;Frederik Barkhof;Philip C. De Witt Hamer;Martin Klein et al., Association between tumor location and neurocognitive functioning using tumor localization maps, , (2019)

[PMID:28246006](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.