



# CD11c Polyclonal Antibody FITC Conjugated

Catalog No: #C91504FITC

Package Size: #C91504FITC 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD11c Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/1163
Conjugates	FITC
Target Name	CD11c
Other Names	CD11c
Accession No.	Uniprot:P20702Gene ID:3687
Calculated MW	126
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

## Published Papers

Wang Yu-Han;Li A-Guo;Wang Hong-Yan;Tu Yong-Sheng et al., Repression of JAK2-STAT1 and PD-L1 by CEP-33779 ameliorates the LPS-induced decline in phagocytic activity of alveolar macrophages and mitigates lung injury in mice, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.