

# phospho-Tau Protein (Thr181) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91641Bio

Package Size: #C91641Bio 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-Tau Protein (Thr181) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Human
Immunogen Type	PK(p-T)PP
Conjugates	Biotin
Target Name	phospho-Tau Protein (Thr181)
Other Names	phospho-Tau Protein (Thr181)
Accession No.	Uniprot:P10636Gene ID:4137
Calculated MW	52/79
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

## Published Papers

Laia Lidon;Cristina Vergara;Isidro Ferrer;Felix Hernandez;Jesus Avila;Jose A. del Rio;Rosalina Gavin et al., Tau Protein as a New Regulator of Cellular Prion Protein Transcription, , (2020)

[PMID:32683652](https://pubmed.ncbi.nlm.nih.gov/32683652/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.