

# Goat Conjugated IgG

Catalog No: #CIG501



Package Size: #CIG501-AF350 100ul #CIG501-AF405 100ul #CIG501-AF488 100ul #CIG501-AF555 100ul #CIG501-AF594 100ul #CIG501-AF647 100ul #CIG501-AF680 100ul #CIG501-AF750 100ul #CIG501-Biotin 100ul  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Conjugated IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,95%
Concentration	20 mg/ml
Storage	4°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.