

Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit

Catalog No: #EK20211



Package Size: #EK20211 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Storage	12 months

Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:7.1pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Published Papers

Zhang, Guojing;Wan, Zilin;Wang, Xiaoyan;Huang, Yunchao;Zhou, Youquan;Chen, Ying;Lei, Yujie;Qiao, Limei et al., In vivo study on the effects of intercellular adhesion operon of Staphylococcus epidermidis on the inflammation associated with bacteria-fungal mixed biofilm, , (2021)

[PMID:34651489](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.