

General Nitrotyrosine (NT) ELISA Kit

Catalog No: #EK20393

Package Size: #EK20393 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	General Nitrotyrosine (NT) ELISA Kit
Applications	ELISA
Species Reactivity	General

Application Details

Detect Range: 7.81-500ng/mL

Sensitivity: 3.5ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Published Papers

Jiao Xueqiao;Liu Moqi;Li Rui;Li Jialu;Wang Lu;Niu Guowei;Wang Liming;Ji Xunming;Lv Chunmei;Guo Xiuhai et al., Helpful to Live Healthier? Intermittent Hypoxic/Ischemic Training Benefits Vascular Homeostasis and Lipid Metabolism with Activating SIRT1 Pathways in Overweight/Obese Individuals, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.