

## APC anti-human CD24

Catalog No: #FC0030



Package Size: #FC0030-1 25Tests #FC0030-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	APC anti-human CD24
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-m15
Isotype	mouse IgG2a, $\kappa$
Applications	FC
Species Reactivity	human
Conjugates	APC
Target Name	CD24
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with APC anti-human CD24 (colored histogram) or APC mouse IgG2a  $\kappa$  (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.