

APC anti-human CD34

Catalog No: #FC0044



Package Size: #FC0044-1 25Tests #FC0044-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


|                    |   |
|--------------------|---|
| Product Name       | APC anti-human CD34   |
| Host Species       | mouse   |
| Clonality          | Monoclonal  |
| Clone No.          | s-4H11  |
| Isotype            | mouse IgG1  |
| Applications       | FC  |
| Species Reactivity | human   |
| Conjugates         | APC   |
| Target Name        | CD34  |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Two mixed TF-1cellsand HEL cellsstained with APC anti-human CD34 (colored histogram) or APC mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.