

APC anti-human CD48

Catalog No: #FC0050



Package Size: #FC0050-1 25Tests #FC0050-2 100Tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	APC anti-human CD48
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-156-4H9
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human
Conjugates	APC
Target Name	CD48
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with APC anti-human CD48 (colored histogram) or APC mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.