

Rabbit IgG Isotype Control FITC Conjugated

Catalog No: #IC001FITC

Package Size: #IC001FITC 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit IgG Isotype Control FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity chromatography purification
Applications	WB, ELISA, FCM, IHC-P, IHC-F, IF(IHC-P), IF(IHC-F), IF(ICC)
Conjugates	conjugated
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol.
Storage	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide. Store at -20°C for 12 months.

Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.