

APC mouse IgG1 α K Isotype Control

Catalog No: #IC002APC



Package Size: #IC002APC-1 25tests #IC002APC-2 100test

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	APC mouse IgG1 α K Isotype Control
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4AIB7.11
Isotype	Mouse IgG1, κ
Purification	Purified by Protein A.
Applications	FC
Formulation	Phosphate-buffered solution, pH 7.2-7.4, containing 0.05% Proclin 300 and 1%FBS
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Vol.per.Test: 10 μ l/Test

Background

This antibody can be used as a mouse IgG1 κ isotype control in flow cytometry and other applications

Note: This product is for in vitro research use only and is not intended for use in humans or animals.