

# PE Mouse IgG1 $\kappa$ Isotype Control

Catalog No: #IC002PE



Package Size: #IC002PE-1 25tests #IC002PE-2 100tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE Mouse IgG1 $\kappa$ Isotype Control
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4AIB7.11
Isotype	Mouse IgG1, $\kappa$
Purification	Purified by Protein A.
Applications	FC IF
Conjugates	conjugated
Concentration	10 $\mu$ g/ml/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA,Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.