

Low Endotoxin Purified rat IgG2b Isotype Control

Catalog No: #IG202



Package Size: #IG202-1 1mg #IG202-2 5mg #IG202-3 25mg #IG202-4 50mg #IG202-5 100mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Low Endotoxin Purified rat IgG2b Isotype Control
Isotype	Rat IgG2b
Purification	>95% by SDS Page
Applications	in vivo
Species Reactivity	ε ^g • iO
Specificity	β 1.0 EU/mg as determined by the LAL method
Formulation	0.01 M η£·ι 'EηO ε ² ζA ² pH 7.2 - 7.4, 150 mM NaCloO ζ θ½½δ½ θ η ½ε ι ² θ • ε
Storage	2-8 εl'οO δΩ• ε- δ'E δ'E ζ οO ι Ωζ δΩ• ε- -80 εl'οO ι• Ωε • ε• • ε • ε »θ •

Application Details

ε ζ— η' ε• «ι • οO ≤ 1.0 EU/mg as determined by the LAL method

Note: This product is for in vitro research use only and is not intended for use in humans or animals.