

# Sheep IgG

Catalog No: #IG301



Package Size: #IG301-1 10mg #IG301-2 1g

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|               |  |
|---------------|--|
| Product Name  | Sheep IgG  |
| Purification  | Affinity chromatography or caprylic acid-ammonium sulfate, $\beta$ 95% |
| Concentration | $\beta$ 10 mg/ml   |
| Storage       | -20°C  |

## Product Description

Source: Normal, healthy, unimmunized sheep serum

Appearance: Liquid

Filtration: 0.2 $\mu$ m

Buffer :10mM sodium Phosphate150mM sodium chloridePBSPH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8 $\Lambda$  C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.