

Goat Anti-Human IgG Secondary Antibody AMCA Conjugated

Catalog No: #L3025

Package Size: #L3025-1 100ul #L3025-2 1.0 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Human IgG Secondary Antibody AMCA Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Human
Specificity	IgG (H+L)
Conjugates	AMCA
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.