

Goat Anti-Mouse IgG Secondary Antibody Rhodamine (TRITC) Conjugated

Catalog No: #L30314

Package Size: #L30314-1 100ul #L30314-2 1.0 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody Rhodamine (TRITC) Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Rhodamine (TRITC)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Published Papers

el at., Clock-modified mesenchymal stromal cells therapy rescues molecular circadian oscillation and age-related bone loss via miR142-3p/Bmal1/YAP signaling axis. In Cell Death Discov. 2022 Mar 12 by Sa Cha 1, Jiangyue Wang, et al.. PMID:35279674, , (2022)

[PMID:35279674](https://pubmed.ncbi.nlm.nih.gov/35279674/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.