

Goat Anti-Rabbit IgG Secondary Antibody AF790 Conjugated

Catalog No: #L30510

Package Size: #L30510-1 100ul #L30510-2 0.5 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody AF790 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rabbit
Specificity	IgG (H+L)
Conjugates	AF647, Alexa 488, FITC
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

Western Blot: 1:20,000-100,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.