

Goat Anti-Rabbit IgG Secondary Antibody AP Conjugated

Catalog No: #L3073



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #L3073-1 100ul #L3073-2 1.0 ml

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody AP Conjugated
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	WB;IHC;ICC;ELISA
Species Reactivity	Rabbit
Specificity	IgG (H+L)
Conjugates	AP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Store at -20 °C/1 year

Application Details

WB 1:1000-1:20000;IHC 1:500-1:1000;ICC 1:500-1:1000;ELISA 1:10000 - 1:50000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.