

Rabbit Anti-Guinea Pig IgG (H&L) Biotin

Catalog No: #L35050



Package Size: #L35050-1 100 ul #L35050-2 500 ul #L35050-3 1 ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Guinea Pig IgG (H&L) Biotin
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Gp
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Guinea Pig IgG
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.