

Goat Anti-Chicken IgY (H&L) FITC

Catalog No: #L35065



Package Size: #L35065-1 100 ul #L35065-2 500 ul #L35065-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Goat Anti-Chicken IgY (H&L) FITC
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Published Papers

el at., SynCAM1 deficiency in the hippocampal parvalbumin interneurons contributes to sevoflurane-induced cognitive impairment in neonatal ratsIn  
CNS Neurosci Ther On2023 Dec 17byMing-Ming Zhao , Han-Wen Gu et al..PMID:38105652, , (2023)  
[PMID:38105652](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.