

Goat Anti-Rabbit IgG (H&L) - AF488

Catalog No: #L35193



Package Size: #L35193-1 100 ul #L35193-2 500 ul #L35193-3 1 ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG (H&L) - AF488
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Rabbit
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Rabbit IgG
Conjugates	AF488, FITC, PE
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Published Papers

Li et al., Targeting STING-mediated pro-inflammatory and pro-fibrotic effects of alveolar macrophages and fibroblasts blunts silicosis caused by silica particlesIn J Hazard MaterOn2023 Sep 15byLiang Ou , Peng Zhang et al..PMID:37379600, , (2023)
[PMID:37379600](https://pubmed.ncbi.nlm.nih.gov/37379600/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.