

# Hairy Cell Leukemia Antibody

Catalog No: #P1056

Package Size: #P1056 0.1m

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Hairy Cell Leukemia Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC687
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

## Application Details

Recommended

working dilution range:1:100 - 1:200

## Images



Hairy Cell Leukemia [P1056] on Lymphoma

## Product Description

Anti-Hairy Cell Leukemia stains various B-cells in the follicular mantle zone and virtually all cases of hairy cell leukemia. It also stains some high grade B-cell lymphomas.

## References

Saati A, et al. Blood. 1989; 74:2476-85.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.