

Recombinant Protein A

Catalog No: #RE006

Package Size: #RE006-1 10mg #RE006-2 1g

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Protein A
Storage	Stability of storage: Recombinant protein A lyophilized should be stored under 2-8°C in sealed container. It is stable within 24 months. When gluing even with PBS (pH7.4) dissolved, it is now with the current use. After dissolved with pure water, it should be stored under -20°C. Stability of transport: The product is stable by blue ice insulation transport.

Product Description

Protein A is a membrane protein of *Staphylococcus aureus* and can specifically bind to antibodies. Protein A affinity chromatography columns have become widely used as affinity columns for purified antibodies from ascites, serum and cell culture isolated or purified extracts of antibodies from different subtypes of mammals or genetically engineered recombinant proteins containing antibody Fc fragments. alkali protein A is our company independent research and development production, with independent intellectual property rights. Recombinant protein A is expressed by *E. coli* and purified by chromatography to give a purity of 95% or more. The purification process does not use antibody column affinity chromatography to avoid the possibility of incorporation of unrelated IgG in the product.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.