

mCherry Mouse Polyclonal Antibody

Catalog No: #T513

Package Size: #T513-1 50ul #T513-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	mCherry Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IF IP
Specificity	The mCherry tag antibody can recognize mCherry and mCherry tag fusion proteins.
Target Name	mCherry
Concentration	1.0mg/ml
Formulation	Mouse IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:5000~1:10000

Images



Western blot analysis of mCherry Recombinant protein, using #T513 diluted at 1) 1:5,000, 2) 1:10,000.

Background

mCherry is a fluorophore (a fluorescent molecule) used in biotechnology as a tracer to follow the flow of fluids, as a marker when tagged to molecules and cells components. mCherry is a monomeric fluorescent construct with peak absorption/emission at 587nm and 610nm, respectively. It is resistant to photobleaching and is stable. mCherry is sometimes preferred to other fluorophores due to its colour, as well as its photostability compared to other monomeric fluorophores.

Published Papers

el at., GABA Cells in the Central Nucleus of the Amygdala Promote Cataplexy. In J Neurosci on 2017 Apr 12 by Matthew B Snow, Jimmy J Fraigne, et

al..PMID: 28209737, , (2017)

PMID:28209737

Note: This product is for in vitro research use only and is not intended for use in humans or animals.