

Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody

Catalog No: #V1017



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1017-1 50ug #V1017-2 100ug


Description

Product Name	Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1649
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM, IF, IHC, SPR, WB
Species Reactivity	Influenza A virus (A/Minnesota/11/2010(H3N2))
Conjugates	Unconjugated
Target Name	Neuraminidase, 3.2.1.18, NA
Accession No.	A0A075ETL7
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-Influenza A virus Neuraminidase/NA Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.