

Anti-Influenza A virus (H1N1) Neuraminidase/NA
Recombinant Antibody

Catalog No: #V1019

Package Size: #V1019-1 50ug #V1019-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Influenza A virus (H1N1) Neuraminidase/NA Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1785
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Influenza A virus (H1N1)
Conjugates	Unconjugated
Target Name	Neuraminidase, EC:3.2.1.18, NA
Accession No.	P03470
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus Neuraminidase/NA Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.