

InVivoMAb Anti-Influenza A virus(strain A/Swine/New Jersey/11/1976 H1N1) HA/Hemagglutinin Antibody

Catalog No: #V1024

Package Size: #V1024-1 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Influenza A virus(strain A/Swine/New Jersey/11/1976 H1N1) HA/Hemagglutinin Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1463
Isotype	IgG1, lambda
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA
Species Reactivity	Influenza A virus (strain A/Swine/New Jersey/11/1976 H1N1)
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	P03455
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Influenza A virus HA/Hemagglutinin Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.