

Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Recombinant Antibody

Catalog No: #V1032

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #V1032-1 50ug #V1032-2 100ug

Description

Product Name	Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1857
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, in vivo protection, in vitro neutralization
Species Reactivity	Influenza A virus H7 subtype
Conjugates	Unconjugated
Target Name	Hemagglutinin (H7 subtype)
Accession No.	P03458, P36346, P26095, P26097, P26101, P26103, P09343, P09344, P26099, P26094, P26094, Q6LEJ4, P26096, P26100
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.