

InVivoMAb Anti-Influenza A
virus(A/Singapore/H2011.447/2011(H3N2))
HA/Hemagglutinin Broadly Neutralizing Antibody

Catalog No: #V1035

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1035-1 100ug


Description

Product Name	InVivoMAb Anti-Influenza A virus(A/Singapore/H2011.447/2011(H3N2)) HA/Hemagglutinin Broadly Neutralizing Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1048
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Influenza A virus (A/Singapore/H2011.447/2011(H3N2)) (R9U684), Influenza A virus (strain A/Victoria/3/1975 H3N2) (P03435), Influenza A virus (A/Victoria/361/2011(H3N2)) (A0A097PF39)
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	R9U684
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



Detects recombinant H3N2 HA in indirect ELISA



SDS-PAGE for InVivoMAb Anti-Influenza A virus
HA/Hemagglutinin Broadly Neutralizing Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.