

Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody

Catalog No: #V1041

Package Size: #V1041-1 50ug #V1041-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1515
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Influenza A virus (A/Vietnam/1194/2004(H5N1))
Conjugates	Unconjugated
Target Name	HA, Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain
Accession No.	Q6DQ34
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.