

Anti-Influenza A virus M/Matrix protein 2 Recombinant Antibody

Catalog No: #V1054

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #V1054-1 50ug #V1054-2 100ug


Description

Product Name	Anti-Influenza A virus M/Matrix protein 2 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1395
Isotype	IgG1-kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IF, WB
Species Reactivity	Influenza A virus (strain A/Puerto Rico/8/1934 H1N1)
Conjugates	Unconjugated
Target Name	Matrix protein 2, Proton channel protein M2, M
Accession No.	P06821
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SEC-HPLC detection for Anti-Influenza A virus M/Matrix protein 2 Antibody .



SDS-PAGE for Anti-Influenza A virus M/Matrix protein 2 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.