

Anti-HPV16 E6/Protein E6 Recombinant Antibody

Catalog No: #V1058



Package Size: #V1058-1 50ug #V1058-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                    |                                                                                                                                                                 |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-HPV16 E6/Protein E6 Recombinant Antibody                                                                                                                   |
| Host Species       | Mouse                                                                                                                                                           |
| Clonality          | Monoclonal                                                                                                                                                      |
| Clone No.          | VS1880                                                                                                                                                          |
| Isotype            | VHH-mFC                                                                                                                                                         |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                             |
| Applications       | ELISA, Neutralization, WB                                                                                                                                       |
| Species Reactivity | Human papillomavirus type 16                                                                                                                                    |
| Conjugates         | Unconjugated                                                                                                                                                    |
| Target Name        | E6, Protein E6                                                                                                                                                  |
| Accession No.      | P03126                                                                                                                                                          |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                |
| Storage            | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks).<br>Store at -20°C 12 months. Store at -80°C long term. |

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-HPV16 E6/Protein E6 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.