

Anti-Human papillomavirus type 16 E6/Protein E6
Recombinant Nanobody

Catalog No: #V1059

Package Size: #V1059-1 50ug #V1059-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Human papillomavirus type 16 E6/Protein E6 Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	VS1871
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA, IF, WB
Species Reactivity	Human papillomavirus type 16
Conjugates	Unconjugated
Target Name	Protein E6, E6
Accession No.	P03126
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-Human papillomavirus type 16 E6/Protein E6 Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.