

Anti-HPV16 E6/Protein E6 Recombinant Antibody

Catalog No: #V1060



Package Size: #V1060-1 50ug #V1060-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV16 E6/Protein E6 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2216
Isotype	VHH-mFc
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Human papillomavirus type 16
Conjugates	Unconjugated
Target Name	E6, Protein E6
Accession No.	P03126
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-HPV16 E6/Protein E6 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.