

Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody

Catalog No: #V1062

Package Size: #V1062-1 50ug #V1062-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1944
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, WB
Species Reactivity	Human papillomavirus type 16
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	P03101
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-HPV16 L1/Major capsid protein L1 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.