

Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody

Catalog No: #V1071



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1071-1 50ug #V1071-2 100ug


Description

Product Name	Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2173
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neut, WB
Species Reactivity	Human papillomavirus type 16 (HPV16)
Conjugates	Unconjugated
Target Name	Major capsid protein L1, L1
Accession No.	P03101
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-HPV16 L1/Major capsid protein L1 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.