

Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody

Catalog No: #V1072



Package Size: #V1072-1 50ug #V1072-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV16 L1/Major capsid protein L1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1973
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Human papillomavirus type 16 (HPV16)
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	P03101
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-HPV16 L1/Major capsid protein L1 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.