

Anti-HPV16 E7/Protein E7 Recombinant Antibody

Catalog No: #V1075

Package Size: #V1075-1 50ug #V1075-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-HPV16 E7/Protein E7 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2135
Isotype	VHH-mFC
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, WB
Species Reactivity	Human papillomavirus type 16
Conjugates	Unconjugated
Target Name	E7, Protein E7
Accession No.	P03129
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.