

Anti-HPV16 L2/Minor capsid protein L2 Recombinant Antibody

Catalog No: #V1084

Package Size: #V1084-1 50ug #V1084-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Anti-HPV16 L2/Minor capsid protein L2 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2268
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	WB
Species Reactivity	Human papillomavirus type (HPV16, HPV18, HPV31, HPV33, HPV35, HPV51, HPV52, HPV58)
Conjugates	Unconjugated
Target Name	Minor capsid protein L2, L2
Accession No.	P03107
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details


Purity:>95% as determined by SDS-PAGE.WB: 1:500-1:2000

Images

Anti-HPV16 L2/Minor capsid protein L2 Antibody (VS2268) binds with HPV16 L2



Detects Human papillomavirus type 16 L2/Minor capsid protein L2 in indirect ELISAs.



SDS-PAGE for Anti-HPV16 L2/Minor capsid protein L2
Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.