

InVivoMAb Anti-HPV16 L2/Minor capsid protein L2 Antibody

Catalog No: #V1085

Package Size: #V1085-1 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-HPV16 L2/Minor capsid protein L2 Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1745
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, WB
Species Reactivity	Human papillomavirus type 16
Conjugates	Unconjugated
Target Name	Minor capsid protein L2, L2
Accession No.	P03107
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



Detects recombinant HPV16 L2/Minor capsid protein L2 in indirect ELISA



SDS-PAGE for InVivoMAb Anti-HPV16 L2/Minor capsid protein L2 .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.