

Anti-HPV18 E7/Protein E7 Recombinant Antibody

Catalog No: #V1087



Package Size: #V1087-1 50ug #V1087-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV18 E7/Protein E7 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2096
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IP, WB, ICC, IHC
Species Reactivity	HPV16, HPV18
Conjugates	Unconjugated
Target Name	E7
Accession No.	P06788
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.

Application Details

Purity:>95% purity as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-HPV18 E7/Protein E7 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.