

Anti-HPV18 E7/Protein E7 Recombinant Antibody

Catalog No: #V1088



Package Size: #V1088-1 50ug #V1088-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV18 E7/Protein E7 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1282
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IHC, WB
Species Reactivity	Human papillomavirus type 18
Conjugates	Unconjugated
Target Name	E7, Protein E7
Accession No.	P06788
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.

Application Details

Purity:>95% purity as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-HPV18 E7/Protein E7 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.