

Anti-HPV33 E7/Protein E7 Recombinant Antibody

Catalog No: #V1096



Package Size: #V1096-1 50ug #V1096-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                    |                                                                                                                                                                                                 |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-HPV33 E7/Protein E7 Recombinant Antibody                                                                                                                                                   |
| Host Species       | Mouse                                                                                                                                                                                           |
| Clonality          | Monoclonal                                                                                                                                                                                      |
| Clone No.          | VS1409                                                                                                                                                                                          |
| Isotype            | IgG2a, kappa                                                                                                                                                                                    |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                                                             |
| Applications       | ELISA, IHC, WB                                                                                                                                                                                  |
| Species Reactivity | Human papillomavirus 33                                                                                                                                                                         |
| Conjugates         | Unconjugated                                                                                                                                                                                    |
| Target Name        | E7, Protein E7                                                                                                                                                                                  |
| Accession No.      | P06429                                                                                                                                                                                          |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                                                |
| Storage            | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C. |

Application Details

Purity:>95% purity as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-HPV33 E7/Protein E7 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.